

GENOMICS RESEARCH CORE

Taqman® SNP Sample Submission**Sample Preparation**

Samples must be submitted in 96 well plates with well H12 empty. Genomic DNA should be at 10ng/ul. If the DNA will be submitted at any other concentration, then this must be noted and discussed with the Genomics Research Core (labor charges will apply). Please [contact](#) the Genomics Research Core for instructions for submission of amplified DNA. Minimum sample volume is 40 μ l for up to 8 assays. Submit an additional 2 ml for each additional assay to be included.

Sample and Assay Submission

Samples and assays can be submitted at any of our drop off locations.

- S534 Scaife Hall
- 3343 Forbes Avenue, 3rd floor
- Room 2.15 Hillman Cancer Center Research Pavilion
- Children's Hospital loading dock
- Room 238, 300 Technology Drive

Please complete a Full Service Realtime PCR sample submission [form](#).

Turn Around Time

Project turn around time will depend on the work queue when samples and assays are delivered as well as size of the project. Please [contact](#) Genomics Research Core at the time of sample submission for estimated date of data delivery.