

GENOMICS RESEARCH CORE

Realtime PCR Full Service Submission

1. Contact the [Genomics Research Core](#) to discuss your project.
2. Obtain assays for all genes.
3. Submit assays and RNA (or biological specimens if Genomics Research Core is to perform RNA extraction) to one of our drop off locations
 - a. Drop off locations
 1. 3343 Forbes Ave Room 303
 2. S534A Scaife Hall
 3. Children's Hospital Loading Dock
 4. Rm 2.15 Hillman Cancer Pavilion
 5. Rm 238 300 Technology Dr.
 - b. Please be sure to complete a [Full Service Taqman Sample Submission Form](#). Hard copies are available in the drop off rooms.
4. You will receive an email when your project is complete. Data will be available from our web accessible folder
 - a. Log in to [MyLogin](#)
 - b. Click on the Other Data link under Get Data
 - c. Click on the link for the current data file.
 - d. Save locally.
5. Genomics Research Core works on a first come first served basis. Therefore the exact turn around for any submission will depend on the amount of work ready and waiting to be performed when your samples arrive. Please consult with Genomics Core staff at the time of project submission regarding the current queue and expected turnaround time.